

LH-TG16E Micro High Speed Centrifuge



Features:

1. Brushless motor, no pollution, free-maintenance.
2. Microprocessor control, LED display which indicates the speed, time, RCF in operation, 10 kinds of brake setting, operate simply.
3. Electric lid lock, super speed, imbalance protection. The centrifuge body is made of high quality steel, stainless steel chamber, safe and reliable.
4. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction.
5. 3 tiers protection steel cover and get the ideal centrifugation result

LH-TG16E Technical Parameter:

Max. Speed	16000rpm	RPM/RCF Convert	Yes
Max. RCF	19040×g	Noise (dB)	≤ 65
Max. Capacity	6×10ml	Temperature	Normal
Time Range	0~99min		

Matched Rotors for LH-TG16E

Order No	Rotor	Max speed (rpm)	Max Volume(ml)	Max RCF (g)
16E-1	Angle rotor	13000	4×8PCR	10410
16E-2	Angle rotor	16000	40×0.2ml	19040
16E-3	Angle Rotor	13000	40×0.5ml	17210
16E-4	Angle Rotor	16000	12×1.5/2ml	17940
16E-5	Angle Rotor	14000	24×1.5/2ml	15470
16E-6	Angle Rotor	13000	10×5ml	11810
16E-7	Angle Rotor	13000	6×10ml	14190
16E-8	Angle Rotor	12000	24 pieces capillary vessel	15800

